


ศรีตรังโกลฟส์®

SRI TRANG GLOVES
SURGICAL

LATEX POWDER-FREE 
SURGICAL GLOVES



Ergonomic design

Conforms to hand contour



Tested against EN 16523-1 and ASTM D6978

for chemical and chemo drugs resistance



NATUTEX™ latex

our soft latex technology proprietary that comes with lower protein content



LATEX POWDER-FREE SURGICAL GLOVES



Material

Natural Rubber Latex

Powder Content

≤ 2.0 mg/glove

Internal Glove Surface

Polymer Coated

Freedom from Holes (Inspection level I)

AQL ≤ 0.65

Product Compliances

ASTM D3577, EN 455,
ISO 10282, TIS 538,
EN ISO 21420, EN ISO 374

Factory Compliances

ISO 13485, ISO 14001,
ISO 9001

Viral Penetration Test

ASTM F1671

PRODUCT INFORMATION

Maximum Protection:

Latex Surgical Gloves offer maximum protection against bloodborne pathogens and viruses, ensuring a safe surgical environment for both surgeons and patients.

Extra Comfort for Prolonged Procedures:

Designed with the comfort of surgeons in mind, these gloves enable long hours of use without compromising comfort or dexterity.

Enhanced Textile Sensitivity:

The special formula compound of Latex enhances textile sensitivity, *allowing surgeons to feel even the most delicate movements*.

Precise Grip:

These gloves provide a precise grip, maximizing surgical tool handling and facilitating precise maneuvers.

Sustainable and Environmentally Friendly:

Latex Surgical Gloves are made from sustainable raw materials and produced using **100% renewable energy**, minimizing their impact on the environment.

Sterile by Gamma Ray:

The gloves undergo *sterilization by gamma ray*, ensuring the highest level of sterility and safety.

Product Dimension

Peel pack (MV)	132(W)x158(L)
Dispenser (Inner box)	137(W)x157(L)x277(H)
Carton (Outer carton)	335(W)x430(L)x287(H)



Container Loading

20ft	40ft	40ft HC
684	1422	1602



Sbk_Mklocal@sritranggroup.com

www.sritranggloves.com