



COMFORT

Latex Powdered Examination Gloves



Reliable protection against **viral penetration**



Micro finger textured
for precise grip



Inner lightly powdered with **USP absorbable dusting powder**

Color



Latex White



Food contact compliance

240 mm length

Ability to reseal

Efficient and easy donning

Single wall thickness (finger)
typical value 4.3 mil (0.11 mm)

Balanced of strength, comfort & flexibility

High elasticity palm reduces hand fatigue

PPE CAT I



Experience the superior comfort of

COMFORT

Latex Powdered Examination Gloves

where hand protection meets unparalleled comfort and sensitivity.

100%
Renewable Energy

Low Carbon Production

Sustainable Sourcing

USAGE
(highlight)



Medical



Hospitality



Beauty



Retail



Food



COMFORT

Latex Powdered Examination Gloves

The ultimate choice for those who prioritize **comfort without compromising on protection.**

These gloves offer a balance between comfort, flexibility, and reliable performance, making them perfect for light tasks that do not require heavy-duty protection.

COMFORT Latex Powdered Examination Gloves are designed to provide a comfortable wearing experience throughout extended periods of use. The latex material offers a soft and supple feel, ensuring optimal comfort and reduced hand fatigue. These gloves are ideal for tasks that require sensitivity and precision, such as light medical procedures, laboratory work, or general tasks where dexterity is essential.

With their thinner profile, **COMFORT** Latex Powdered Examination Gloves offer enhanced tactile sensitivity and a second-skin feel. This allows for precise movements and increased control, making them a preferred choice for delicate tasks where a high level of tactile feedback is required.

Furthermore, *these gloves are lightly powdered with USP absorbable dusting powder*, ensuring easy donning and smooth glove application. The powdered interior reduces friction, allowing for effortless glove changes and enhancing overall user convenience.

COMFORT Latex Powdered Examination Gloves are suitable for various industries and applications where a lighter level of protection is needed. They are well-suited for light medical tasks, dental procedures, laboratory work, beauty salons, and other environments that require sensitivity, comfort, and adequate hand protection.

| | | | |
|--|--|--|---------------------------|
| <p>Length</p> <p>Average thickness at palm (Single wall)</p> <p>Tensile Strength</p> <p>Force at Break</p> <p>Elongation</p> | | | <p>E-Document QR code</p> |
| | <p>XS-S : Min 220mm, M-XL : Min 230mm</p> | <p>Median \geq 240mm</p> | |
| | <p>Min 0,05mm</p> | <p>0,09mm</p> | |
| | <p>BEFORE AGING AFTER AGING</p> <p>18MPa 14MPa</p> | <p>BEFORE AGING AFTER AGING</p> <p>20MPa 19MPa</p> | |
| | <p>BEFORE AGING AFTER AGING</p> <p>650% 500%</p> | <p>BEFORE AGING AFTER AGING</p> <p>710% 670%</p> | |

glove_global@sritranggroup.com

www.sritranggloves.com

KEY FEATURES:



Comfort-Oriented Design:

Prioritizes comfort and reduces hand fatigue, providing a comfortable wearing experience throughout extended use.



Enhanced Tactile Sensitivity:

Thinner profile offers increased tactile sensitivity and a second-skin feel, allowing for precise movements and improved control.



Lightly Powdered for Easy Donning:

Inner lightly powdered with **USP absorbable dusting powder**, ensuring smooth and convenient glove application.



Versatile Application:

Suitable for light medical procedures, dental tasks, laboratory work, beauty salons, and other light-duty applications that require sensitivity and comfort.



Reliable Performance:

Provides adequate hand protection for tasks that do not require heavy-duty protection, offering peace of mind during various light-duty tasks.



Sustainable and Environmentally Friendly:

Made from sustainable raw materials and produced using **100% renewable energy**, reducing environmental impact

