



PREMIUM

Latex Powdered Examination Gloves



Reliable protection against **viral penetration**



Micro finger textured
for precise grip



Inner lightly powdered with **USP absorbable dusting powder**



Color

Latex White

Single wall thickness (finger) **typical value**
5.5 mil (0.14 mm)

Excellent experience **donning**

High elasticity
palm reduces hand fatigue

240 mm length

Food contact compliance

Ability to reseal

Superior strength, comfort & flexibility

MDR CLASS 1



PPE CAT III
CE 2777



PPE CAT I
CE



Experience the trusted quality of

PREMIUM

Latex Powdered Examination Gloves

where comfort, protection, sustainability, and hygiene standards align for optimal performance.

100%
Renewable Energy

Low Carbon
Production

Sustainable
Sourcing

USAGE
(highlight)



Medical



EMS



Industry



Food



Dental



PREMIUM

Latex Powdered Examination Gloves

The classic glove series that has been trusted and widely used for over a decade.

These gloves offer exceptional thickness, comfort, and flexibility, with an elongation capability over 700%. They provide a snug fit that allows for natural hand movements and enhanced dexterity, ensuring optimal performance in various tasks.

PREMIUM Latex Powdered Examination Gloves are the superb choice when it comes to protection against viral penetration. Their reliable thickness and high-quality latex material create a robust barrier that effectively safeguards against potential hazards. These gloves are designed to meet the stringent requirements of hygiene standards, making them ideal for medical, dental, food handling, and general hand protection applications.

Ease of donning is a key feature of **PREMIUM Latex Powdered Examination Gloves**. They are *lightly powdered with USP absorbable dusting powder*, ensuring smooth and hassle-free glove application. This feature enhances user convenience, making these gloves a preferred choice for tasks that demand efficiency and quick glove changes.

In addition to their protective properties and user convenience, **PREMIUM Latex Powdered Examination Gloves** prioritize sustainability. They are produced using *100% renewable energy*, minimizing their environmental impact and contributing to a greener future. By choosing these gloves, you not only prioritize the safety and well-being of yourself and others but also support eco-friendly practices.



Comply to



E-Document QR code

Length

XS-S : Min 220mm, M-XL : Min 230mm

Median ≥ 240mm

Median ≥ 240mm

Average thickness at palm (Single wall)

Min 0.05mm

-

0.11mm

Tensile Strength

18MPa | 14MPa

-

28MPa | 26MPa

Force at Break

-

6.0N | 6.0N

8.0N | 7.5N

Elongation

650% | 500%


-

780% | 750%

✉ glove_global@sritranggroup.com

🌐 www.sritranggloves.com

KEY FEATURES:




Enhanced Grip:

Finger-textured surface provides excellent grip, allowing for better control and precision in tasks that require a secure hold.



Superb Protection:

Reliable defense against viral penetration, ensuring a high level of protection in various environments.



Lightly Powdered for Easy Donning:

Inner lightly powdered with *USP absorbable dusting powder*, ensuring smooth and convenient glove application.



Classic Gloves Series:

Trusted and widely used for over a decade, known for their exceptional *thickness, comfort, and flexibility*.



Sustainable and Environmentally Friendly:

Made from sustainable raw materials and produced using *100% renewable energy*, reducing environmental impact.