



COMFORT

Nitrile Powder-Free Examination Gloves




Micro texture at fingertip
precise grip



High comfort & flexibility



Great protection
against viral penetration and chemical permeation



Color

Ocean Blue

Violet Blue

White

Efficient and easy donning

240 mm length

Food contact compliance



Single wall thickness (finger) typical value 3.5 mil (0.09 mm)

MDR CLASS 1



PPE CAT III
CE 2777



PPE CAT I
CE



Experience
the difference of

COMFORT

Nitrile Powder-Free Examination Gloves

and enjoy enhanced comfort and performance throughout your tasks.

100%
Renewable Energy



Low Carbon
Production

USAGE
(highlight)



Medical



Dental



Laboratory



Food



Retail



COMFORT

Nitrile Powder-Free Examination Gloves

The perfect choice for individuals seeking reliable protection combined with exceptional **comfort and reduced hand fatigue.**

These gloves are designed to provide a comfortable and reliable solution for a wide range of applications, making them ideal for medical professionals and anyone working in environments that require protective equipment and safe food handling practices.

COMFORT Nitrile Examination Gloves excel in various fields within the medical industry and beyond. Whether you're performing medical examinations, handling protective equipment, or ensuring food safety, these gloves offer reliable performance and peace of mind.

One of the key features of **COMFORT Nitrile Examination Gloves** is their exceptional comfort and skin-like feel. The soft and comfortable nitrile material ensures a pleasant wearing experience, reducing discomfort and fatigue during extended use. This is particularly beneficial for tasks that require prolonged glove use, minimizing hand strain and promoting better productivity.

In addition to providing outstanding comfort, **COMFORT Nitrile Examination Gloves** offer good chemical resistance. The gloves are designed to withstand exposure to a range of chemicals commonly encountered in medical settings, ensuring reliable protection for the wearers.

KEY FEATURES:



Reliable Protection and Comfort:

Combines reliable protection with exceptional comfort, reducing hand fatigue and **promoting a comfortable wearing experience.**



Versatile Applications:

Suitable for a wide range of applications, including medical examinations, protective equipment handling, and safe food handling.



Skin-Like Feel:

The soft and comfortable nitrile material provides a skin-like feel, enhancing wearer comfort and **reducing discomfort during extended use.**



Reduced Hand Fatigue:

Minimizes hand strain and fatigue, allowing for more comfortable and prolonged glove use.



Good Chemical Resistance:

Withstands exposure to a range of chemicals commonly encountered in medical environments, ensuring reliable protection for both the wearer and the patient.



Produced Using 100% Renewable Energy:

Made with a commitment to sustainability, these gloves are produced using **100% renewable energy**, reducing their environmental impact.

Comply to			
Length	XS-S : Min 220mm, M-XL : Min 230mm	Median ≥ 240mm	Median ≥ 240mm
Average thickness at palm (Single wall)	Min 0.05mm	-	0.07mm
Tensile Strength	14MPa 14MPa	-	40MPa 44MPa
Force at Break	-	6.0N 6.0N	6.8N 7.5N
Elongation	500% 400%	-	610% 550%

Scan for more information / Scan for product details / Scan for sustainability

E-Document QR code

glove_global@sritranggroup.com

www.sritranggloves.com

