



# SHIELD PLUS Nitrile Powder Free Examination Gloves



**Tested for use against chemotherapy drugs**  
(acc. to ASTM D6978)



**Tested against EN 16523-1 and ASTM F739**  
for chemical resistance



**No permeation against Fentanyl up to 240 minutes**  
(acc. to ASTM D6978)



**Color**

Blue

Single wall thickness (finger) **typical value 5.1 mil (0.13 mm)**

Efficient and **easy donning**

Powerful hand protection against **multiple hazards**

 **Food contact compliance**

**300 mm length**

**Extended cuff** for extended protection

**Micro textured finger** precise grip

**Great chemical resistance**

**MDR CLASS 1**



**PPE CAT III**  
CE 2777



**PPE CAT I** CE



**300 mm LONG CUFF**  
**EXTRA PROTECTION**

With their extra-long cuff, exceptional durability, and outstanding chemical resistance, these gloves offer unparalleled defense and versatility.

**100%**  
Renewable Energy



**Low Carbon**  
Production

**USAGE**  
(highlight)



Hospital



Dental



Laboratory



Janitorial



Industry



Food



# SHIELD PLUS Nitrile Powder Free Examination Gloves

The ultimate choice for those seeking extra protection and versatility. These gloves are specially designed with a **300mm extra-long cuff**, providing unparalleled coverage and ensuring full protection against chemicals and liquids entering through the wrist area. Whether you're working with hazardous substances or wearing a PPE suit, these gloves are the perfect complement, offering comprehensive protection from top to bottom.

SHIELD PLUS Nitrile Examination Gloves strike an impressive balance between thickness and flexibility, making them a reliable choice for various tasks. They offer exceptional durability and protection, surpassing the capabilities of thin nitrile gloves. From medium-duty work to more demanding tasks, these gloves provide reliable defense against cuts, abrasions, and punctures, instilling confidence and peace of mind.

One of the standout features of SHIELD PLUS Nitrile Examination Gloves is their outstanding chemical resistance. They act as a formidable barrier against a wide range of chemicals, safeguarding your hands during chemical-intensive tasks. These gloves have been tested for **chemotherapy drug protection** and have demonstrated **no permeation against fentanyl for up to 240 minutes**. This level of protection makes them an ideal choice for healthcare professionals, laboratory technicians, and those working with hazardous substances.

In addition to their exceptional protection, SHIELD PLUS Nitrile Examination Gloves are highly versatile. They are suitable for both industrial applications and food handling, offering a reliable solution for various work environments. From manufacturing plants to food processing facilities, these gloves provide the necessary protection and performance required to carry out tasks safely and efficiently.



glove\_global@sritrangroup.com

www.sritrangloves.com

KEY FEATURES:




### Extra Protection with 300mm Long Cuff:

Ensures full protection from chemicals and liquids entering through the wrist area, providing comprehensive coverage.



### Balanced Thickness and Flexibility:

Strike the perfect balance between thickness and flexibility, offering optimal durability, protection, and comfort.



### Chemical Resistant and Tested for Chemotherapy Drug Protection:

Formidable barrier against a wide range of chemicals, including protection against chemotherapy drugs and no permeation against fentanyl for up to 240 minutes.



### Versatile Performance:

Suitable for both industrial and food handling applications, catering to a wide range of work environments.



### Suitable for Medical Use:

Compliant with rigorous quality standards, making them an excellent choice for healthcare professionals.

